

ST3GAL3 antibody



Catalog Number: 115623

Product name

ST3GAL3 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ST3GAL3 Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue. Observed molecular weight by Western-blot: 44-55 kDa

Preparation

This antibody was obtained by immunization of ST3GAL3 recombinant protein (Accession Number: NM_006279). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:115623(ST3GAL3 Antibody) at dilution of 1:300