

# ST7 antibody



Catalog Number: 115631

## Product name

ST7 antibody

Recommended Dilution:

WB: 1:1000-1:10000

## Specificity

Human, Mouse; other species not tested.

IP: 1:500-1:5000

IHC: 1:20-1:200

## Antibody description

ST7 Rabbit Polyclonal antibody. Positive IHC detected in human liver tissue, human heart tissue, human lymphoma tissue, human pancreas tissue. Positive IP detected in PC-3 cells. Positive WB detected in mouse brain tissue, HeLa cells, human brain tissue, Jurkat cells, K-562 cells, PC-3 cells. Observed molecular weight by Western-blot: 63 kDa

## Preparation

This antibody was obtained by immunization of ST7 recombinant protein (Accession Number: BC030954). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type


Rabbit IgG

## Applications


ELISA, IHC, IP, WB

## Dilutions

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:115631(ST7 Antibody) at dilution of 1:2000




IP Result of anti-ST7 (IP:Catalog No:115631, 4ug; Detection:Catalog No:115631 1:1500) with PC-3 cells lysate 2640ug.


# ST7 antibody

Catalog Number: 115631

---



Immunohistochemical of paraffin-embedded human liver using Catalog No:115631(ST7 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human liver using Catalog No:115631(ST7 antibody) at dilution of 1:50 (under 40x lens)