

# ST8SIA2-Specific antibody



Catalog Number: 115634

## Product name

ST8SIA2-Specific antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ST8SIA2-Specific Rabbit Polyclonal antibody.  
Positive WB detected in mouse heart tissue.  
Positive IHC detected in human heart tissue.  
Observed molecular weight by Western-blot: 42kd

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_006011).  
Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol  
pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200


## Validations



mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:115634(ST8SIA2-Specific antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human heart using Catalog No:115634(ST8SIA2-Specific antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human heart using Catalog No:115634(ST8SIA2-Specific antibody) at dilution of 1:50 (under 40x lens)