

# STAC2 antibody

Catalog Number: 115637

## Product name

STAC2 antibody

## Specificity

Human; other species not tested.

## Antibody description

STAC2 Rabbit Polyclonal antibody. Positive WB detected in SH-SY5Y cells. Positive IHC detected in human prostate hyperplasia tissue. Observed molecular weight by Western-blot: 42-45 kDa

## Preparation

This antibody was obtained by immunization of STAC2 recombinant protein (Accession Number: NM\_001351360). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

## Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:115637(STAC2 Antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded human prostate hyperplasia slide using Catalog No:115637(STAC2 Antibody) at dilution of 1:50