

SPR antibody

Catalog Number: 115641

Product name

SPR antibody

IP: 1:200-1:2000

IF: 1:10-1:100


Specificity

Human, Mouse, Rat; other species not tested.

Validations

Antibody description

SPR Rabbit Polyclonal antibody. Positive IF detected in Hela cells. Positive IP detected in mouse liver tissue. Positive WB detected in HepG2 cells, human liver tissue, K-562 cells, rat heart tissue, rat liver tissue. Observed molecular weight by Western-blot: 28 kDa




Preparation

This antibody was obtained by immunization of SPR recombinant protein (Accession Number: NM_003124). Purification method: Antigen affinity purified.

HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:115641 (SPR antibody) at dilution of 1:500

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Immunofluorescent analysis of Hela cells, using SPR antibody Catalog No:115641 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).

Ig Type

Rabbit IgG


Applications

ELISA, WB, IF, IP

Dilutions

Recommended Dilution:

WB: 1:500-1:5000



IP Result of anti-SPR (IP:Catalog No:115641, 4ug; Detection:Catalog No:115641 1:500) with mouse

SPR antibody



Catalog Number: 115641

liver tissue lysate 4000ug.