

SPTBN2 antibody



Catalog Number: 115654

Product name

SPTBN2 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

SPTBN2 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, mouse kidney tissue. Observed molecular weight by Western-blot: 270kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_006946). Purification method: Antigen affinity purified.

Formulation

0.1M NaHCO₃, 0.1M glycine, 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:115654(SPTBN2 antibody) at dilution of 1:500