

# SREBF1 antibody

Catalog Number: 115669

## Product name

SREBF1 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

SREBF1 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, L02 cells. Positive IP detected in L02 cells. Observed molecular weight by Western-blot: 125kd

## Preparation

This antibody was obtained by immunization of SREBF1 recombinant protein (Accession Number: NM\_001005291). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:115669(SREBF1 antibody) at dilution of 1:500



IP Result of anti-SREBF1 (IP:Catalog No:115669, 4ug; Detection:Catalog No:115669 1:600) with L02 cells lysate 1500ug.