

# STARD3NL antibody



Catalog Number: 115679

## Product name

STARD3NL antibody

## Immunogen

[Human STARD3NL Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

STARD3NL Rabbit Polyclonal antibody. Positive IF detected in HEK-293 cells. Positive IHC detected in human breast cancer tissue, human gliomas tissue. Positive IP detected in HEK-293 cells. Positive WB detected in HEK-293 cells, COLO 320 cells. Observed molecular weight by Western-blot: 27 kDa

## Preparation

This antibody was obtained by immunization of STARD3NL recombinant protein (Accession Number: XM\_017012694). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC, IP, IF

## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IP: 1:200-1:2000


IHC: 1:20-1:200

IF: 1:20-1:200

## Validations




HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:115679(STARD3NL antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:115679(STARD3NL antibody) at dilution of 1:100 (under 10x lens)


# STARD3NL antibody

Catalog Number: 115679




Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:115679(STARD3NL antibody) at dilution of

1:100 (under 40x lens)



IP Result of anti-STARD3NL (IP:Catalog No:115679, 4ug; Detection:Catalog No:115679 1:500) with HEK-293 cells lysate 3200ug.



Immunofluorescent analysis of (10% Formaldehyde) fixed HEK-293 cells using Catalog No:115679(STARD3NL Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)