

STAT4 antibody



Catalog Number: 115689

Product name

Polyclonal

STAT4 antibody

Ig Type

Specificity

Rabbit IgG

human, mouse, rat; other species not tested.

Applications

Antibody description

ELISA, WB

STAT4 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HeLa cells, PC-3 cells. Observed molecular weight by Western-blot: 86kd

Dilutions

Recommended Dilution:

Preparation


WB: 1:200-1:2000

This antibody was obtained by immunization of STAT4 recombinant protein (Accession Number: XM_006712719). Purification method: Antigen affinity purified.

Validations

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:115689(STAT4 antibody) at dilution of 1:300