

STAT6 antibody



Catalog Number: 115696

Product name

STAT6 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

STAT6 Rabbit Polyclonal antibody. Positive IHC detected in human spleen. Positive WB detected in NIH/3T3 cells, MCF7 cells. Observed molecular weight by Western-blot: 110kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: XM_047429475). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with Catalog No:115696(STAT6 antibody) at dilution of 1:1500



Immunohistochemical of paraffin-embedded human spleen using Catalog No:115696(STAT6 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human spleen using Catalog No:115696(STAT6 antibody) at dilution of 1:50 (under 40x lens)