

STAU1 antibody



Catalog Number: 115699

Product name

STAU1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

STAU1 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in human brain tissue, K-562 cells, mouse brain tissue. Observed molecular weight by Western-blot: 55 kDa

Preparation

This antibody was obtained by immunization of STAU1 recombinant protein (Accession Number: NM_001322930). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:115699(STAU1 antibody) at dilution of 1:300



IP Result of anti-STAU1 (IP:Catalog No:115699, 4ug; Detection:Catalog No:115699 1:500) with mouse brain tissue lysate 3440ug.