

STC1 antibody



Catalog Number: 115702

Product name

STC1 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

STC1 Rabbit Polyclonal antibody. Positive WB detected in mouse small intestine tissue. Observed molecular weight by Western-blot: 31kd

Preparation

This antibody was obtained by immunization of STC1 recombinant protein (Accession Number: NM_003155). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse small intestine tissue were subjected to SDS PAGE followed by western blot with Catalog No:115702(STC1 antibody) at dilution of 1:500