

# STIP1 antibody



Catalog Number: 115709

## Product name

STIP1 antibody

## Immunogen

[Human STIP1 Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

STIP1 Rabbit Polyclonal antibody. Positive IP detected in Jurkat cells. Positive WB detected in Jurkat cells. Positive IF detected in MCF-7 cells. Observed molecular weight by Western-blot: 61 kDa

## Preparation

This antibody was obtained by immunization of STIP1 recombinant protein (Accession Number: NM\_006819). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF, IP

## Dilutions


Recommended Dilution:

WB: 1:200-1:2000


IP: 1:200-1:2000

IF: 1:10-1:100

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:115709(STIP1 antibody) at dilution of 1:800



Immunofluorescent analysis of MCF-7 cells, using STIP1 antibody Catalog No:115709 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit

# STIP1 antibody

Catalog Number: 115709

---

IP Result of anti-STIP1 (IP:Catalog No:115709, 4ug; Detection:Catalog No:115709 1:500) with Jurkat cells lysate 2400ug.

