

# STK3 antibody

Catalog Number: 115715

## Product name

STK3 antibody

IHC: 1:20-1:200

IF: 1:20-1:200


## Specificity

Human, Mouse, Rat; other species not tested.

## Validations

## Antibody description

STK3 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, A375 cells. Positive IF detected in HeLa cells. Positive IHC detected in human ovary tumor tissue. Observed molecular weight by Western-blot: 63 kDa




## Preparation

This antibody was obtained by immunization of STK3 recombinant protein (Accession Number: NM\_006281). Purification method: Antigen affinity purified.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:115715(STK3 antibody) at dilution of 1:300

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human ovary tumor using Catalog No:115715(STK3 antibody) at dilution of 1:50 (under 10x lens)

## Ig Type

Rabbit IgG


## Applications

ELISA, WB, IHC, IF

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000



Immunofluorescent analysis of HeLa cells, using

# STK3 antibody



Catalog Number: 115715

---

STK3 antibody Catalog No:115715 at 1:50 dilution and FITC-labeled donkey anti-rabbit IgG (green).  
Blue pseudocolor = DAPI (fluorescent DNA dye).