

STOML3 antibody



Catalog Number: 115732

Product name

STOML3 antibody

Immunogen

[Human STOML3 Recombinant protein \(His tag\)](#)

Specificity

human, mouse, rat; other species not tested.

Antibody description

STOML3 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, MCF7 cells, SH-SY5Y cells. Observed molecular weight by Western-blot: 30kd

Preparation

This antibody was obtained by immunization of STOML3 recombinant protein (Accession Number: NM_145286). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:500-1:5000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:115732(STOML3 antibody) at dilution of 1:500