

# STRAP antibody

Catalog Number: 115735

## Product name

STRAP antibody

WB: 1:1000-1:10000

IP: 1:1000-1:10000

## Specificity


Human, Mouse, Rat; other species not tested.

IHC: 1:20-1:200

## Antibody description

STRAP Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, HeLa cells, HepG2 cells. Positive IP detected in HeLa cells. Positive IHC detected in human ovary tissue, human lung tissue, human placenta tissue, human skin tissue, human testis tissue. Observed molecular weight by Western-blot: 39 kDa

## Validations



## Preparation

This antibody was obtained by immunization of STRAP recombinant protein (Accession Number: NM\_007178). Purification method: Antigen affinity purified.

mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:115735(STRAP antibody) at dilution of 1:2000

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type


Rabbit IgG

## Applications


ELISA, WB, IHC, IP

## Dilutions

Recommended Dilution:



IP Result of anti-STRAP (IP:Catalog No:115735, 4ug; Detection:Catalog No:115735 1:2000) with HeLa cells lysate 2000ug.




# STRAP antibody

Catalog Number: 115735

---

Immunohistochemical of paraffin-embedded human ovary using Catalog No:115735(STRAP antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human ovary using Catalog No:115735(STRAP antibody) at dilution of 1:100 (under 40x lens)