

# STRBP antibody

Catalog Number: 115736

## Product name

STRBP antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

STRBP Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells, HeLa cells, mouse testis tissue. Positive IP detected in PC-3 cells. Observed molecular weight by Western-blot: 74 kDa

## Preparation

This antibody was obtained by immunization of STRBP recombinant protein (Accession Number: XM\_017014894). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IP: 1:200-1:2000

## Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:115736(STRBP antibody) at dilution of 1:500



IP Result of anti-STRBP (IP:Catalog No:115736, 4ug; Detection:Catalog No:115736 1:500) with PC-3 cells lysate 4000ug.