

SYN2 antibody



Catalog Number: 115764

Product name

SYN2 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

SYN2 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 63 kDa

Preparation

This antibody was obtained by immunization of SYN2 recombinant protein (Accession Number: NM_133625). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:115764(SYN2 Antibody) at dilution of 1:2000