

Synaptogyrin 1 antibody



Catalog Number: 115765

Product name

Synaptogyrin 1 antibody

Specificity

Human,Mouse,Rat; other species not tested.

Antibody description

Synaptogyrin 1 Rabbit Polyclonal antibody.
Positive WB detected in mouse brain tissue.
Observed molecular weight by Western-blot: 28-30 kDa

Preparation

This antibody was obtained by immunization of Synaptogyrin 1 recombinant protein (Accession Number: NM_145731). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:115765(SYNGR1 antibody) at dilution of 1:100