

Synaptotagmin-3 antibody



Catalog Number: 115775

Product name

Synaptotagmin-3 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Synaptotagmin-3 Rabbit Polyclonal antibody.
Positive WB detected in mouse thymus tissue.
Observed molecular weight by Western-blot: 63kd

Preparation

This antibody was obtained by immunization of Synaptotagmin-3 recombinant protein (Accession Number: NM_001160328). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse thymus tissue were subjected to SDS PAGE followed by western blot with Catalog No:115775(SYT3 antibody) at dilution of 1:300