

# Synaptotagmin-4 antibody



Catalog Number: 115777

## Product name

Synaptotagmin-4 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

Synaptotagmin-4 Rabbit Polyclonal antibody.  
Positive WB detected in SH-SY5Y cells, A431 cells,  
mouse brain tissue. Observed molecular weight by  
Western-blot: 48-50 kDa

## Preparation

This antibody was obtained by immunization of  
Synaptotagmin-4 recombinant protein (Accession  
Number: NM\_020783). Purification method:  
Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol  
pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



SH-SY5Y cells were subjected to SDS PAGE  
followed by western blot with Catalog  
No:115777(SYT4 antibody) at dilution of 1:500