

SYNC antibody

Catalog Number: 115779

Product name

SYNC antibody

Immunogen

[Human SYNC Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

SYNC Rabbit Polyclonal antibody. Positive IP detected in mouse heart tissue. Positive WB detected in mouse heart tissue. Positive IHC detected in human skeletal muscle tissue, human heart tissue. Observed molecular weight by Western-blot: 64-68, 55 kDa

Preparation

This antibody was obtained by immunization of SYNC recombinant protein (Accession Number: BC119700). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


IP: 1:200-1:2000

IHC: 1:20-1:200

Validations




mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:115779(SYNC Antibody) at dilution of 1:600




IP Result of anti-SYNC (IP:Catalog No:115779, 4ug; Detection:Catalog No:115779 1:800) with mouse heart tissue lysate 4000ug.

SYNC antibody

Catalog Number: 115779



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue slide using Catalog No:115779(SYNC Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue slide using Catalog No:115779(SYNC Antibody) at dilution of 1:50 (under 40x lens)