

SYNJ2BP antibody



Catalog Number: 115785

Product name

SYNJ2BP antibody

Immunogen

[Human ARIP2; SYNJ2BP Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

SYNJ2BP Rabbit Polyclonal antibody. Positive IF detected in MCF-7 cells. Positive IHC detected in human stomach cancer tissue. Positive IP detected in mouse lung tissue. Positive WB detected in human heart tissue, A549 cells, human kidney tissue, mouse lung tissue. Observed molecular weight by Western-blot: 16kd

Preparation

This antibody was obtained by immunization of SYNJ2BP recombinant protein (Accession Number: NM_018373). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IF, IP

Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IP: 1:200-1:2000


IHC: 1:20-1:200

IF: 1:10-1:100

Validations




human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:115785(ARIP2; SYNJ2BP antibody) at dilution of 1:500




Immunohistochemical of paraffin-embedded human stomach cancer using Catalog No:115785(ARIP2; SYNJ2BP antibody) at dilution of 1:50 (under 40x lens)

SYNJ2BP antibody

Catalog Number: 115785



Immunofluorescent analysis of MCF-7 cells, using SYNJ2BP antibody Catalog No:115785 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit



IgG (red).

IP Result of anti-ARIP2; SYNJ2BP (IP:Catalog No:115785, 3ug; Detection:Catalog No:115785 1:500) with mouse lung tissue lysate 4000ug.