

Syntaxin 17 antibody

Catalog Number: 115794

Product name

Syntaxin 17 antibody

WB: 1:500-1:5000

IP: 1:500-1:5000

Specificity

Human, Mouse, Rat; other species not tested.

IHC: 1:20-1:200

IF: 1:20-1:200


Antibody description

Syntaxin 17 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive IHC detected in human liver tissue. Positive IP detected in HepG2 cells. Positive WB detected in L02 cells, HepG2 cells, human testis tissue, mouse testis tissue. Observed molecular weight by Western-blot: 33 kDa

Validations

Preparation


This antibody was obtained by immunization of Syntaxin 17 recombinant protein (Accession Number: XM_011518821). Purification method: Antigen affinity purified.



Immunohistochemical of paraffin-embedded human liver using Catalog No:115794(STX17 antibody) at dilution of 1:50 (under 10x lens)

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Immunohistochemical of paraffin-embedded human liver using Catalog No:115794(STX17 antibody) at dilution of 1:50 (under 40x lens)

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC, WB, IP, IF

Dilutions


Recommended Dilution:

Syntaxin 17 antibody


Catalog Number: 115794



IP Result of anti-STX17 (IP:Catalog No:115794, 4ug; Detection:Catalog No:115794 1:1000) with HepG2 cells lysate 3600ug.



L02 cells were subjected to SDS PAGE followed by western blot with Catalog No:115794(STX17 Antibody) at dilution of 1:1000



Immunofluorescent analysis of HepG2 cells using Catalog No:115794(STX17 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)