

# Syntaxin 4 antibody

Catalog Number: 115798

## Product name

Syntaxin 4 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

Syntaxin 4 Rabbit Polyclonal antibody. Positive IP detected in mouse kidney tissue. Positive WB detected in Jurkat cells, HEK-293 cells, mouse kidney tissue, rat kidney tissue. Observed molecular weight by Western-blot: 35-37 kDa

## Preparation

This antibody was obtained by immunization of Syntaxin 4 recombinant protein (Accession Number: NM\_004604). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:115798 (STX4 antibody) at dilution of 1:1500



IP Result of anti-STX4 (IP:Catalog No:115798, 4ug; Detection:Catalog No:115798 1:1500) with mouse kidney tissue lysate 4000ug.