

SULF2 antibody

Catalog Number: 115824

Product name

SULF2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SULF2 Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer tissue, human ovary tissue. Positive WB detected in mouse liver tissue. Observed molecular weight by Western-blot: 130-150 kDa

Preparation

This antibody was obtained by immunization of SULF2 recombinant protein (Accession Number: BC020962). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:115824(SULF2 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:115824(SULF2 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human breast cancer using Catalog

SULF2 antibody



Catalog Number: 115824

No:115824(SULF2 antibody) at dilution of 1:50 (under 40x lens)