

# TALDO1 antibody



Catalog Number: 115844

---

## Product name

TALDO1 antibody

## Specificity

Human,Mouse,Rat; other species not tested.

## Antibody description

TALDO1 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, COLO 320 cells, human placenta tissue. Observed molecular weight by Western-blot: 38kd

## Preparation

This antibody was obtained by immunization of TALDO1 recombinant protein (Accession Number: BC009680). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:115844(TALDO1 antibody) at dilution of 1:500