

# TAOK3 antibody

Catalog Number: 115849

## Product name

TAOK3 antibody

IP: 1:200-1:2000

IHC: 1:20-1:200

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

TAOK3 Rabbit Polyclonal antibody. Positive IHC detected in human liver tissue. Positive IP detected in HepG2 cells. Positive WB detected in human liver tissue, HeLa cells, HepG2 cells, human brain tissue. Observed molecular weight by Western-blot: 100-105 kDa

## Preparation

This antibody was obtained by immunization of TAOK3 recombinant protein (Accession Number: XM\_047428948). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB, IHC, IP

## Dilutions


Recommended Dilution:

WB: 1:200-1:2000


## Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:115849(TAOK3 antibody) at dilution of 1:800



Immunohistochemical of paraffin-embedded human liver using Catalog No:115849(TAOK3 antibody) at dilution of 1:50 (under 10x lens)




Immunohistochemical of paraffin-embedded human liver using Catalog No:115849(TAOK3

# TAOK3 antibody

Catalog Number: 115849

---

antibody) at dilution of 1:50 (under 40x lens)



IP Result of anti-TAOK3 (IP:Catalog No:115849, 4ug; Detection:Catalog No:115849 1:700) with HepG2 cells lysate 3600ug.