

TAU antibody

Catalog Number: 115857

Product name

TAU antibody

WB: 1:500-1:5000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TAU Rabbit Polyclonal antibody. Positive WB detected in SH-SY5Y cells, HeLa cells, HT-1080 cells, mouse brain tissue, rat brain tissue. Positive IHC detected in human brain tissue, human gliomas tissue, human lung cancer tissue. Observed molecular weight by Western-blot: 68 kDa, 81 kDa, 100 kDa

Preparation

This antibody was obtained by immunization of TAU recombinant protein (Accession Number: NM_016841). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:115857(TAU antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human brain tissue slide using Catalog No:115857(TAU Antibody) at dilution of 1:200 (under 10x lens)



Immunohistochemistry of paraffin-embedded

TAU antibody



Catalog Number: 115857

human brain tissue slide using Catalog

No:115857(TAU Antibody) at dilution of 1:200
(under 40x lens)