

# TAX1BP1 antibody



Catalog Number: 115858

## Product name

TAX1BP1 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

TAX1BP1 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, HepG2 cells, mouse kidney tissue, mouse liver tissue, mouse lung tissue, mouse testis tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 68 kDa, 85 kDa

## Preparation

This antibody was obtained by immunization of TAX1BP1 recombinant protein (Accession Number: NM\_006024). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:115858(TAX1BP1 antibody) at dilution of



Immunofluorescent analysis of HepG2 cells using Catalog No:115858(TAX1BP1 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)