

TBC1D1 antibody

Catalog Number: 115861

Product name

TBC1D1 antibody

IP: 1:500-1:5000

IHC: 1:20-1:200

Specificity

Human, Mouse; other species not tested.

Antibody description

TBC1D1 Rabbit Polyclonal antibody. Positive IHC detected in human testis tissue, human kidney tissue. Positive WB detected in Jurkat cells, HeLa cells. Positive IP detected in Jurkat cells. Observed molecular weight by Western-blot: 150-165 kDa

Preparation

This antibody was obtained by immunization of TBC1D1 recombinant protein (Accession Number: XM_047449890). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:500-1:5000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:115861(TBC1D1 antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human testis slide using Catalog No:115861(TBC1D1 Antibody) at dilution of 1:50



Immunohistochemistry of paraffin-embedded human testis slide using Catalog No:115861(TBC1D1 Antibody) at dilution of 1:50

TBC1D1 antibody

Catalog Number: 115861

IP Result of anti-TBC1D1 (IP:Catalog No:115861, 4ug; Detection:Catalog No:115861 1:1000) with Jurkat cells lysate 2400ug.

