

TBL1XR1 antibody



Catalog Number: 115876

Product name

TBL1XR1 antibody

Polyclonal

Ig Type

Rabbit IgG

Specificity

Human, Mouse; other species not tested.

Applications

ELISA, WB

Antibody description

TBL1XR1 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, MCF-7 cells. Observed molecular weight by Western-blot: 55-65 kDa

Dilutions


Recommended Dilution:

WB: 1:200-1:1000

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_024665). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:115876(TBL1XR1 Antibody) at dilution of 1:300

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality