

# TBP antibody

Catalog Number: 115878

## Product name

TBP antibody

## Immunogen

[Human TBP Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse; other species not tested.

## Antibody description

TBP Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, COLO 320 cells, HeLa cells, mouse colon tissue, SKOV-3 cells. Positive IF detected in MCF-7 cells. Observed molecular weight by Western-blot: 46kd

## Preparation

This antibody was obtained by immunization of TBP recombinant protein (Accession Number: NM\_003194). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:115878(TBP antibody) at dilution of 1:1000



Immunofluorescent analysis of MCF-7 cells, using TBP antibody Catalog No:115878 at 1:50 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).