

# TBP antibody

Catalog Number: 115879

## Product name

TBP antibody

## Immunogen

[Human TBP Recombinant protein \(His tag\)](#)

## Specificity

Human, Mouse; other species not tested.

## Antibody description

TBP Rabbit Polyclonal antibody. Positive IP detected in COLO 320 cells. Positive WB detected in COLO 320 cells. Observed molecular weight by Western-blot: 43 kDa

## Preparation

This antibody was obtained by immunization of TBP recombinant protein (Accession Number: NM\_003194). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP

## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

## Validations



COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:117339(TBP antibody) at dilution of 1:1000



IP Result of anti-TBP (IP:Catalog No:117339, 4ug; Detection:Catalog No:117339 1:500) with COLO 320 cells lysate 400ug.