

# TBRG1 antibody



Catalog Number: 115883

## Product name

TBRG1 antibody

## Immunogen

[Human TBRG1 Recombinant protein \(His tag\)](#)

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

TBRG1 Rabbit Polyclonal antibody. Positive WB detected in BxPC-3 cells. Observed molecular weight by Western-blot: 45kd

## Preparation

This antibody was obtained by immunization of TBRG1 recombinant protein (Accession Number: BC032312). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:2000

## Validations



BxPC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:115883(TBRG1 antibody) at dilution of 1:600