

TBX15 antibody



Catalog Number: 115886

Product name

TBX15 antibody

Immunogen

[Human TBX15 Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

TBX15 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 55 kDa

Preparation

This antibody was obtained by immunization of TBX15 recombinant protein (Accession Number: NM_152380). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:115886(TBX15 Antibody) at dilution of 1:600