

# TBX2 antibody

Catalog Number: 115889

## Product name

TBX2 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

TBX2 Rabbit Polyclonal antibody. Positive WB detected in mouse lung tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 74 kDa

## Preparation

This antibody was obtained by immunization of TBX2 recombinant protein (Accession Number: BC052566). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:20-1:200

## Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:115889(TBX2 Antibody) at dilution of 1:300



Immunofluorescent analysis of HepG2 cells using Catalog No:115889(TBX2 Antibody) at dilution of 1:50 and Rhodamine-Goat anti-Rabbit IgG