

TBX21,T-bet antibody

Catalog Number: 115890

Product name

TBX21,T-bet antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

TBX21,T-bet Rabbit Polyclonal antibody. Positive WB detected in mouse thymus tissue. Positive IP detected in mouse thymus tissue. Observed molecular weight by Western-blot: 60-62 kDa

Preparation

This antibody was obtained by immunization of TBX21,T-bet recombinant protein (Accession Number: NM_013351). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:1000

Validations



mouse thymus tissue were subjected to SDS PAGE followed by western blot with Catalog No:115890(TBX21,T-bet antibody) at dilution of



1:300

IP Result of anti-TBX21,T-bet (IP:Catalog No:115890, 4ug; Detection:Catalog No:115890 1:300) with mouse thymus tissue lysate 2800ug.