

TBX22 antibody

Catalog Number: 115891

Product name

TBX22 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TBX22 Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells, mouse testis tissue. Observed molecular weight by Western-blot: 45 kDa and 55 kDa

Preparation

This antibody was obtained by immunization of TBX22 recombinant protein (Accession Number: NM_016954). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:115891(TBX22 antibody) at dilution of 1:500