

# TBX3 antibody

Catalog Number: 115892

## Product name

TBX3 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

TBX3 Rabbit Polyclonal antibody. Positive IP detected in human placenta tissue. Positive WB detected in mouse placenta tissue, HepG2 cells. Observed molecular weight by Western-blot: 80 kDa

## Preparation

This antibody was obtained by immunization of TBX3 recombinant protein (Accession Number: NM\_016569). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

## Validations



mouse placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:115892(TBX3 antibody) at dilution of 1:500



IP Result of anti-TBX3 (IP:Catalog No:115892, 4ug; Detection:Catalog No:115892 1:500) with human placenta tissue lysate 1520ug.