

TBX5 antibody

Catalog Number: 115893

Product name

TBX5 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TBX5 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in HeLa cells, HepG2 cells, MCF7 cells. Observed molecular weight by Western-blot: 58kd

Preparation

This antibody was obtained by immunization of TBX5 recombinant protein (Accession Number: NM_181486). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:115893(TBX5 antibody) at dilution of 1:300



Immunofluorescent analysis of HepG2 cells using Catalog No:115893(TBX5 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)