

# TCEAL7 antibody



Catalog Number: 115899

## Product name

TCEAL7 antibody

## Specificity

human; other species not tested.

## Antibody description

TCEAL7 Rabbit Polyclonal antibody. Positive WB detected in human liver tissue, A549 cells. Positive IHC detected in human ovary tumor. Observed molecular weight by Western-blot: 12 kDa

## Preparation

This antibody was obtained by immunization of TCEAL7 recombinant protein (Accession Number: NM\_152278). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC

## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:115899(TCEAL7 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human ovary tumor using Catalog No:115899(TCEAL7 antibody) at dilution of 1:50 (under 40x lens)