

TCEAL7 antibody



Catalog Number: 115899

Product name

TCEAL7 antibody

Specificity

human; other species not tested.

Antibody description

TCEAL7 Rabbit Polyclonal antibody. Positive WB detected in human liver tissue, A549 cells. Positive IHC detected in human ovary tumor. Observed molecular weight by Western-blot: 12 kDa

Preparation

This antibody was obtained by immunization of TCEAL7 recombinant protein (Accession Number: NM_152278). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:115899(TCEAL7 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human ovary tumor using Catalog No:115899(TCEAL7 antibody) at dilution of 1:50 (under 40x lens)