

# TCF12/HEB antibody



Catalog Number: 115904

## Product name

TCF12/HEB antibody

IHC: 1:20-1:200

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

TCF12/HEB Rabbit Polyclonal antibody. Positive WB detected in mouse lung tissue, Raji cells. Positive IHC detected in human cervical cancer tissue. Observed molecular weight by Western-blot: 73kd

## Preparation

This antibody was obtained by immunization of TCF12/HEB recombinant protein (Accession Number: NM\_207037). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB, IHC

## Dilutions


Recommended Dilution:

WB: 1:200-1:1000


## Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:115904(TCF12 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human cervical cancer using Catalog No:115904(TCF12 antibody) at dilution of 1:50 (under 40x lens)



Immunohistochemical of paraffin-embedded human cervical cancer using Catalog

# TCF12/HEB antibody



Catalog Number: 115904

---

No:115904(TCF12 antibody) at dilution of 1:50 (under 10x lens)