

TCN1 antibody

Catalog Number: 115912

Product name

TCN1 antibody

Specificity

Human; other species not tested.

Antibody description

TCN1 Rabbit Polyclonal antibody. Positive IHC detected in human testis tissue, human tonsillitis tissue. Positive WB detected in human saliva tissue. Observed molecular weight by Western-blot: 65-70 kDa

Preparation

This antibody was obtained by immunization of TCN1 recombinant protein (Accession Number: BC018632). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IHC: 1:20-1:200


Validations



human saliva tissue were subjected to SDS PAGE followed by western blot with Catalog No:115912(TCN1 Antibody) at dilution of 1:2000



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:115912(TCN1 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:115912(TCN1 Antibody) at dilution of 1:50

TCN1 antibody



Catalog Number: 115912

(under 40x lens)