

TCN2 antibody

Catalog Number: 115913

Product name

TCN2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TCN2 Rabbit Polyclonal antibody. Positive IHC detected in human lung cancer tissue. Positive WB detected in A375 cells, mouse small intestine tissue, rat liver tissue, Transfected HEK-293 cells. Observed molecular weight by Western-blot: 48kd

Preparation

This antibody was obtained by immunization of TCN2 recombinant protein (Accession Number: NM_001184726). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



A375 cells were subjected to SDS PAGE followed by western blot with Catalog No:115913(TCN2 antibody) at dilution of 1:100



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:115913(TCN2 antibody) at dilution of 1:50 (under 10x lens)