

# TDG antibody



Catalog Number: 115922

## Product name

TDG antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

TDG Rabbit Polyclonal antibody. Positive IP detected in U-937 cells. Positive WB detected in human brain tissue, HeLa cells, mouse colon tissue, mouse thymus tissue, Raji cells. Observed molecular weight by Western-blot: 60, 86 kDa

## Preparation

This antibody was obtained by immunization of TDG recombinant protein (Accession Number: NM\_003211). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:1000

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:115922(TDG antibody) at dilution of 1:400



IP Result of anti-TDG (IP:Catalog No:115922, 3ug; Detection:Catalog No:115922 1:400) with U-937 cells lysate 3000ug.