

# TERF2IP antibody



Catalog Number: 115948

## Product name

TERF2IP antibody

Recommended Dilution:

WB: 1:1000-1:10000

## Specificity

Human, Mouse, Rat; other species not tested.

IP: 1:200-1:2000


IHC: 1:20-1:200

## Antibody description

TERF2IP Rabbit Polyclonal antibody. Positive IHC detected in human tonsillitis tissue. Positive IF detected in HeLa cells. Positive FC detected in HeLa cells. Positive IP detected in HeLa cells. Positive WB detected in HEK-293 cells, A431 cells, HeLa cells, HT-1080 cells, human placenta tissue, K-562 cells. Observed molecular weight by Western-blot: 60kd

IF: 1:10-1:100

## Validations



## Preparation

This antibody was obtained by immunization of TERF2IP recombinant protein (Accession Number: NM\_018975). Purification method: Antigen affinity purified.

HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:115948(TERF2IP antibody) at dilution of 1:800

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal


## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF, IHC, IP, FC

## Dilutions




Immunofluorescent analysis of HeLa cells, using TERF2IP antibody Catalog No:115948 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit

# TERF2IP antibody




Catalog Number: 115948




IgG (red).

IP Result of anti-TERF2IP (IP:Catalog No:115948, 3ug; Detection:Catalog No:115948 1:800) with HeLa cells lysate 2720ug.



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:115948(TERF2IP Antibody) at dilution of 1:50 (under 40x lens)



1X10<sup>6</sup> HeLa cells were stained with 0.2ug TERF2IP antibody (Catalog No:115948, red) and control antibody (blue). Fixed with 90% MeOH blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1000.