

TACC1 antibody

Catalog Number: 115967

Product name

TACC1 antibody

Specificity

Human; other species not tested.

Antibody description

TACC1 Rabbit Polyclonal antibody. Positive IF detected in HeLa cells. Positive WB detected in MCF7 cells, HeLa cells. Observed molecular weight by Western-blot: 110-130kd

Preparation

This antibody was obtained by immunization of TACC1 recombinant protein (Accession Number: NM_001352786). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:115967(TACC1 antibody) at dilution of 1:1000



Immunofluorescent analysis of HeLa cells using Catalog No:115967(TACC1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)