

TACR2 antibody



Catalog Number: 115971

Product name

TACR2 antibody

Specificity

Human; other species not tested.

Antibody description

TACR2 Rabbit Polyclonal antibody. Positive FC detected in PC-3 cells. Positive WB detected in HeLa cells, PC-3 cells. Observed molecular weight by Western-blot: 44 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_001057).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, FC

Dilutions


Recommended Dilution:

WB: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:115971(TACR2 Antibody) at dilution of 1:600



1X10⁶ PC-3 cells were stained with 0.2ug TACR2 antibody (Catalog No:115971, red) and control antibody (blue). Fixed with 4% PFA blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1500.